













































Symbol Glossary

Symbol	Symbol Ref. No.	Symbol Title	Symbol Description	Standard
	5.1.1	Manufacturer	Indicates the medical device manufacturer, as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC.	ISO 15223-1:2016 Medical Devices – Symbols to Be Used with Medical Device Labels, Labeling, and Information to be Supplied
	5.1.2	Authorized Representative in the European Community	Indicates the Authorized Representative in the European Community.	ISO 15223-1:2016 Medical Devices – Symbols to Be Used with Medical Device Labels, Labeling, and Information to be Supplied
	5.1.4	Use-by Date	Indicates the date after which the medical device is not to be used.	ISO 15223-1:2016 Medical Devices – Symbols to Be Used with Medical Device Labels, Labeling, and Information to be Supplied
	5.1.5	Batch Code	Indicates the manufacturer's batch code so that the batch or lot can be identified NOTE: Synonyms for "batch code" are "lot number" and "batch number"	ISO 15223-1:2016 Medical Devices – Symbols to Be Used with Medical Device Labels, Labeling, and Information to be Supplied
	5.1.6	Catalogue Number	Indicates the manufacturer's catalog number so that the medical device can be identified NOTE: Synonyms for "catalogue number" are "reference number" and "reorder number"	ISO 15223-1:2016 Medical Devices – Symbols to Be Used with Medical Device Labels, Labeling, and Information to be Supplied
	5.2.3	Sterilized Using Ethylene Oxide	Indicates a medical device that has been sterilized using ethylene oxide.	ISO 15223-1:2016 Medical Devices – Symbols to Be Used with Medical Device Labels, Labeling, and Information to be Supplied
	5.2.4	Sterilized Using Irradiation	Indicates a medical device that has been sterilized using irradiation.	ISO 15223-1:2016 Medical Devices – Symbols to Be Used with Medical Device Labels, Labeling, and Information to be Supplied
	5.2.6	Do Not Resterilize	Indicates a medical device that is not to be resterilized.	ISO 15223-1:2016 Medical Devices – Symbols to Be Used with Medical Device Labels, Labeling, and Information to be Supplied
	5.2.7	Non-Sterile	Indicates a medical device that has not been subjected to a sterilization process.	ISO 15223-1:2016 Medical Devices – Symbols to Be Used with Medical Device Labels, Labeling, and Information to be Supplied
	5.2.8	Do Not Use if Package is Damaged	Indicates a medical device that should not be used if the package has been damaged or opened.	ISO 15223-1:2016 Medical Devices – Symbols to Be Used with Medical Device Labels, Labeling, and Information to be Supplied
	5.4.2	Do Not Reuse	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.	ISO 15223-1:2016 Medical Devices – Symbols to Be Used with Medical Device Labels, Labeling, and Information to be Supplied
	5.4.3	Consult Instructions for Use	Indicates the need for the user to consult the instructions for use.	ISO 15223-1:2016 Medical Devices – Symbols to Be Used with Medical Device Labels, Labeling, and Information to be Supplied

Symbol	Symbol Ref. No.	Symbol Title	Symbol Description	Standard
	5.4.4	Caution	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.	ISO 15223-1:2016 Medical Devices – Symbols to Be Used with Medical Device Labels, Labeling, and Information to be Supplied
	3.1.11	MR Conditional	An item with demonstrated safety in the MR environment within defined conditions.	ASTM F2503-13 Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment
	2725	Contains or presence of	On medical devices: to indicate that the equipment contains the identified product or substance.	ISO 7000:2019 Graphical symbols for use on equipment -- Registered symbols

Symbol	Symbol Description	Additional Information
	Rx only. Federal (USA) law restricts this device to sale by or on the order of a physician. An alternative to the text "Rx Only."	21 CFR 801.109 Prescription devices
	CE Marking. The product conforms to European Medical Directive 93/42/EEC and meets applicable health, safety and environmental requirements. If the mark is accompanied by a number, conformity is verified by the indicated notified body.	Medical Device Directive 93/42/EEC
	Indicates unit quantity	Symbol used by A&E Medical Corporation
	A&E Single Wire Sternum Suture (Single Suture)	Symbol used by A&E Medical Corporation
	A&E DoubleWire™ Sternum Suture (Double Suture)	Symbol used by A&E Medical Corporation
	1/2 Circle Needle	Symbol used by A&E Medical Corporation
	3/8 Circle Needle (Pacing Wire)	Symbol used by A&E Medical Corporation
	3/8 Circle Needle (Sternum Suture)	Symbol used by A&E Medical Corporation
	Reverse Cutting Needle Tip	Symbol used by A&E Medical Corporation
	Conventional Cutting Needle Tip	Symbol used by A&E Medical Corporation
	Taper Cutting Needle Tip	Symbol used by A&E Medical Corporation
	Diamond Point Needle Tip	Symbol used by A&E Medical Corporation
	Ball End Needle Tip	Symbol used by A&E Medical Corporation
	Taper Point Needle Tip	Symbol used by A&E Medical Corporation
	Open Pouch in this Direction	Symbol used by A&E Medical Corporation

	Indicates standard style equipment connector	Symbol used by A&E Medical Corporation
	Indicates universal pin equipment connector	Symbol used by A&E Medical Corporation
	Do not connect to electrical outlets	Symbol used by A&E Medical Corporation
B	Blue color	Symbol used by A&E Medical Corporation
W	White color	Symbol used by A&E Medical Corporation
O	Orange color	Symbol used by A&E Medical Corporation
LO	Light orange color	Symbol used by A&E Medical Corporation
	Epicardial MYO/WIRE® (M-20)	Symbol used by A&E Medical Corporation
	Epicardial Loop MYO/WIRE® (M-21)	Symbol used by A&E Medical Corporation
	Epicardial Tapered Loop MYO/WIRE® (M-22)	Symbol used by A&E Medical Corporation
	Snap-off Keith Needle	Symbol used by A&E Medical Corporation
	Straight Needle	Symbol used by A&E Medical Corporation
	Intramyocardial MYO/WIRE® (M-24)	Symbol used by A&E Medical Corporation
	Intramyocardial MYO/WIRE® (M-25)	Symbol used by A&E Medical Corporation
	Intramyocardial MYO/WIRE® (M-26)	Symbol used by A&E Medical Corporation
	Indifferent Lead (M-27)	Symbol used by A&E Medical Corporation
	Indifferent Lead (M-27) no curved needle	Symbol used by A&E Medical Corporation
	Number of needles in pouch	Symbol used by A&E Medical Corporation

5206 Asbury Road, PO Box 758
Farmingdale, NJ 07727 USA
Tel: 732-938-2266
www.aemedical.com